



DB131: Ovalbumin (FL12)+HRP

Background:

Ovalbumin is a glycoprotein with molecular weight of 45 kDa. The molecule consists of a polypeptide with up to two phosphate groups per mole and a side chain of mannose and glucosamine residues.

Origin:

Ovalbumin (FL12)+HRP is provided as an affinity purified rabbit polyclonal that has been raised against highly purified chicken egg albumin .

Product Details:

200 µg/0.2 ml rabbit polyclonal IgG conjugated to HRP in 0.2 ml PBS containing 0.01% thimerosal.

Specificity:

Ovalbumin (FL12)+HRP DB131 forms a single precipitin arc with purified chicken egg white ovalbumin when assayed by immunoelectrophoresis. It is suitable for all ELISA and blotting procedures. Recommended starting dilution for Western Blotting is 1:2000.

Storage:

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.