

# DB092: T7 (E11)

## **Background:**

T7 is an epitope tag composed of an 11residue peptide, MASMTGGQQMG, derived from the T7 major capsid protein of the T7 phage.

# **Origin:**

T7 (E11) is provided as an affinity purified rabbit polyclonal antibody, raised against a peptide derived from the T7 major capsid protein of the T7 phage.

#### **Product Details:**

Each vial contains 200  $\mu$ g/ml of affinity purified rabbit IgG, T7 (E11) DB092, in 1 ml PBS containing 0.1 % sodium azide and 0.2% gelatin.

## **Specificity:**

T7 (E11) is recommended for use by western blotting and immunoprecipitation. The suggested western blotting starting dilution is at 1:400.

#### **Storage:**

Store this product at 4° C, do not freeze. The product is stable for one year from the date of shipment.